

Instructions for measuring the wheelchair to use HURT-e

1. Wheelchair identification

The first information is the name of **the manufacturer, the type of wheelchair and the width of the seat** (all catalog values or values from the label placed on the wheelchair). Next, measure the internal distance between the main tubes of the wheelchair (measure with a tape measure).



2. Photo documentation of the wheelchair

To get an idea of the construction of the wheelchair, we need to send photos of the wheelchair from the **side, front, top and bottom view**. These views must be perpendicular to the side of the photo being taken and should always cover the whole wheelchair.



3. Details and dimensions

Other views are already photographed in detail. This is a view of the **space under the seat + detail of the brake attachment** and seat construction - the entire tube parallel to the main wheels. If there are free circular holes in the structure, **measure their diameter and spacing** (measure with a caliper). It will be possible to use them for fastening with screws.



4. Tube diameters

The second mounting option is to place the sleeves on **parallel main tubes** or on a **cross tube** under the seat. Here it is necessary to measure the diameters of the tubes (measure with a caliper). If the design uses oval tubes, measure both dimensions perpendicular to each other. Take a picture of the space between the seat and the main tube + the space from the end of the seat to the bend of the main tube.



5. Front wheels attachment

The third variant is the attachment to the structure of the **front swivel wheels**, either with screws or with sleeves. Take a picture of the detail of the front wheel mount from the outside and inside view. Measure the diameter of the holes / diameter of the screws + their spacing (measure with a caliper).



Please return the required images and dimensions for review to info@hurt-e.cz. According to the provided data, we will prepare a mounting kit for your wheelchair. If necessary, we will request additional data. The assumption is that the assembly is carried out by a skilled craftsman who has the basic tools at his disposal.

Wheelchair measure form

Customer

Name: _____

Contact: _____

Wheelchair identification

Manufacturer: _____ Model: _____

Seat width (label data): _____ mm / Width between main tubes (measured): _____ mm

Details and dimensions

Main tube

Diameter: _____ mm / Oval dimensions: _____ mm, _____ mm

Holes in the tube

Diameter: _____ mm Spacing: _____ mm

Cross tube

Diameter: _____ mm / Oval dimensions: _____ mm, _____ mm

Front swivel wheels

Diameter: _____ mm Spacing: _____ mm

HURT-e color

Black

Silver

White

Yellow "banana"

Red

Orange "Kia"

Blue

Other _____

Technical requirements

Main brake with arrest left or right

Additional battery

Accelerator on the left

One brake lever for both brakes

Notes
