

HURT-e

Wheelchair Power Attachment



Main advantages of HURT-e

- One universal wheelchair mobility device.
- Very light construction.
- Perfectly space saving = small dimensions (in the car trunk, below the bed, ...).
- Easy removable from the wheelchair by quick attachment.
- Great power and range.
- Option to connect a second battery.
- Turn the wheelchair in place by turning the handlebar by 90 °.
- Choice of several levels of power.
- Possibility to purchase with the trolley.

Parameters

- Wheel size - 10 inches.
- Motor - 48 V / 500 W.
- Battery 48 V / 8,8 Ah.
- Charging time - 6 hours.
- Maximum range - 25 km.
- Weight - 12 kg.
- Maximum speed - up to 40 km/h under various conditions.
- Maximum load - depending on wheelchair load.
- Dimensions - 90/50/30 cm (height of the device/handlebar width/wheel length and handle).
- Disc brake.
- Display data: speed, distance, battery capacity...
- Programmable user parameters.



[hurt-e.CZ](http://hurt-e.cz)

HURT s. r. o.
Bambousek 664
281 26 Týnec nad Labem
CZECH REPUBLIC
Tel.: +420 608 977 274
E-mail: info@hurt-e.cz