



BASIC

Basic version with a range of up to 30 km

PRICE: 2790 €



DUO

Double range up to 60 km

PRICE: 3349 €

It is only possible to make preparations for the later purchase of a battery

129 €



TRIO

Triple range up to 90 km

PRICE: 3590 €

Main advantages of HURT-e

- One device for wheelchair user's mobility.
- We are the manufacturers of this device ourselves.
- Very light and quality construction.
- Perfectly storage = small dimensions (in the car trunk, below the bed, behind the wardrobe ...).
- Easy removal from the wheelchair thanks to quick release attachment.
- Great power and range.
- Option to connect a second battery for greater range.
- Turn the wheelchair on the spot by turning the handlebars by 90°.
- Choice of three levels of power.
- For fixed and folding frame wheelchairs (folding retained).
- Possibility to purchase with a wheelchair.
- Included: installation of attachment to the wheelchair, choice of basic colours, charger, casing, rear light for wheelchair.

Parameters

- Wheel size - 10 inches.
- Motor - 48 V / 500 W.
- Battery 48 V (BASIC / DUO / TRIO) 8,8 / 17,6 / 23,8 Ah.
- Charging time (BASIC / DUO / TRIO) 6 / 12 / 18 h.
- Maximum range 30 / 60 / 90 km (depends on many variables).
- Weight (BASIC / DUO / TRIO) 12 / 15 / 17 kg.
- Maximum speed – 25 km/h according to legislation.
- Design speed up to 40 km/h. Setup at your own responsibility with limited warranty.
- Maximum load - depending on wheelchair load.
- Dimensions – 90 / 50 / 30 cm (height of the device / handlebar width / length with attachment).
- Disc brake with two independent brakes.
- Brake lever lock.
- Display data: speed, distance, battery capacity...
- Reverse gear.
- Programmable user parameters.
- Possibility of further individual adjustments.



HURT s. r. o.
Bambousek 664
281 26 Týnec nad Labem
CZECH REPUBLIC

hurt-e.CZ

Tel.: +420 777 327 910
E-mail: info@hurt-e.cz

Accessories

Color on request

79 €



Specify your color according to the RAL swatch.

QUADRU control

809 €



Control the accelerator by tilting the handlebars towards, brake by moving away from yourself. Suitable for quadriplegics.

CROSS tire

12 € / 36 €
primary production/spare part



Easier passage over moderate terrain.

Hydraulic brake

209 €



Great brake effect in combination with a higher quality disc.

Multi-purpose bag

35 €



Easy installation. We can send by mail.

Basket for handlebars

39 €



The metal basket can be easily removed thanks to the simple CLICK system. Easy installation. We can send by mail.

Adjustable fork

31 €



For quadriplegic control. Adjustable.

Left-hand accelerator

129 €



Left-handed throttle control.

Detailed information including photos and videos on hurt-e.cz.

We work with partners



Czech manufacturer of wheelchairs



Special driving school for the disabled