

# HURT-e

# Wheelchair Power Attachment



## Main advantages of HURT-e

- One universal wheelchair mobility device.
- Very light and quality construction.
- Perfectly space saving = small dimensions (in the car trunk, below the bed, ...).
- Easy removable from the wheelchair by quick release attachment.
- Great power and range.
- Option to connect a second battery for greater range.
- Turn the wheelchair in place by turning the handlebar by 90°.
- Choice of three levels of power.
- Possibility to purchase with a wheelchair.

**Price: 2599 EUR** (including our 21 % VAT)

Included: installation of attachment to the wheelchair, choice of basic colours, casing.

## Modernization for model year 2021

- Reverse gear.
- Brake boost - Shimano main brake.
- Selection of the left / right brake lever for the main braking effect.
- Selection of left / right brake lever with lock for additional braking effect.
- Anodized surface finish of the connecting bracket.

## Optional packages

- Second battery for double range.
- Left-hand gas.
- Possibility of only one brake lever with effect on both brakes.
- Color on request.
- Possibility of further individual adjustments.



**hurt-e.CZ**

HURT s. r. o.  
Bambousek 664  
281 26 Týnec nad Labem  
Czech Republic  
Tel.: +420 608 977 274  
E-mail: info@hurt-e.cz



## Parameters

- Wheel size - 10 inches.
- Motor - 48 V / 500 W.
- Battery 48 V / 8,8 Ah.
- Charging time - 6 hours.
- Maximum range - 30 km (depends on many variables).
- Weight - 12 kg.
- Maximum speed – 25 km/h according to legislation.
- Design speed up to 40 km/h. Adjust at your own risk.
- Maximum load - depending on wheelchair load.
- Dimensions - 90/50/30 cm (height of the device/handlebar width/length with attachment).
- Disc brake.
- Display data: speed, distance, battery capacity...

